

Implant System	Model	Size	SmartPeg	Article no.
STRAUMANN DENTAL IMPLANT SYSTEM	Bone Level Narrow CrossFit	3.3	Type 53	100484
	Bone Level Regular CrossFit	4.1, 4.8	Type 54	100485
	Extra Oral		Type 1	100353
	Ortho/Palatal		Type 9	100368
	Roxolid Standard, Standard Plus, TE	Regular Neck 3.3	Type 4	100350
	Roxolid Bone Level Narrow CrossFit	3.3	Type 41	100463
	Standard, Standard Plus, TE	Narrow Neck	Type 17	100390
		Regular Neck	Type 4	100350
		Wide Neck	Type 4	100350
	SWEDEN & MARTINA	Kohno - new platform	3.8	Type 32
Kohno - old platform		3.8	Type 21	100402
Kohno - old + new platform		4.25, 5.0	Type 7	100380
Premium - old platform		3.3, 3.8	Type 21	100402
Premium - old + new platform		4.25, 5.0	Type 7	100380
Global			Type 48	100479
SYBRON	Endopore External	4.1, 5.0	Type 1	100353
	Entegra External	3.25, 4.0	Type 1	100353
	Pitt-Easy	4.0, 4.9	Type 24	100412
	Sybron PRO XRT	3.5, 4.1	Type 6	100378
TBR	External Hexagon, Internal Hexagon		Type 1	100353
	TEKKA	In-Kone, T-Quest	Type 49	100480
	Krestal, Progress		Type 6	100378
	TwinKon		Type 17	100390
THOMMEN MEDICAL	SPI System	platform 3.5	Type 22	100404
		platform 4.0	Type 21	100402
		platform 4.5 - 6.0	Type 15	100386
TIDAL SPIRAL IMPLANT SYST.		3.75, 4.2, 5.0, 6.0	Type 27	100431
U-IMPL	Internal Hexagon		Type 32	100440
	External Hexagon		Type 1	100353
WIELAND DENTAL IMPLANTS		3.5, 4.3, 5.0	Type 45	100472
ZIMMER DENTAL	AdVent	4.5	Type 35	100446
		5.7	Type 34	100444
	Screw-Vent (also Tapered)	3.5	Type 32	100440
		4.5	Type 27	100431
		5.7	Type 33	100442
	Spline	3.25 - 5.0	Type 30	100436
	SwissPlus (also Tapered)	platform 3.8	Type 31	100438
		platform 4.8	Type 4	100350

Abutment System	Model	Size	SmartPeg	Article no.
COCHLEAR BAHA	BA 210		Type 55	100487
	BA 300		Type 55	100487
NEOSS	Access Abutment		Type 52	100483
NOBEL BIOCARE	Brånemark System MUA	NP, RP	Type A3	100363
THOMMEN MEDICAL	SPI SYSTEM Bar Abutment		Type 36	100447
	VARIOmulti Abutment		Type 25	100415

Please note that this document is updated continuously. The latest version can always be found at www.osstell.com.

SmartPeg Reference List

valid from March 2011

Implant System	Model	Size	SmartPeg	Article no.	
3i	Osseotite/Nanotite External Hex MicroMiniplant	3.25	Type 1	100353	
	Osseotite/Nanotite External Hex	3.75, 4.0, 5.0, 6.0	Type 1	100353	
	Osseotite/Nanotite External Hex	XP 3/4, 4/5, 5/6	Type 1	100353	
	Osseotite/Nanotite Certain MicroMiniplant	3.25	Type 45	100472	
	Osseotite/Nanotite Certain Prevail	4/3	Type 45	100472	
	Osseotite/Nanotite Certain	4.0, 5.0, 6.0	Type 15	100386	
	Osseotite/Nanotite Certain Prevail	5/4, 4/5/4, 5/6/5	Type 15	100386	
	Osseotite Certain	XP 4/5, 5/6	Type 15	100386	
	A.B. DENTAL DEVICES	Screw Type, Conical Groovy	3.75	Type 32	100440
			4.5	Type 27	100431
ADIN DENTAL IMPLANT SYSTEMS	One		Type 17	100390	
	Swell		Type 49	100480	
	Touareg, Touareg-S		Type 49	100480	
ALPHA BIO IMPLANT	ATI, SPI		Type 26	100425	
ALFA-GATE	Bioactive		Type 35	100446	
ALPHATECH		3.4 - 7.0	Type 46	100477	
ANTHOGYR	Anthofit OI Tapered, OI Straight		Type 18	100398	
	Anthofit HE Straight		Type 1	100353	
	Axiom	3.4-5.2	Type 38	100455	
	Ossfit		Type 4	100350	
ASTRA TECH	Osseospeed Yellow	3.0 S	Type 50	100481	
	Osseospeed Aqua	3.5 S, 4.0 S	Type 38	100455	
	Osseospeed Lilac	4.5, 5.0, 5.0 S	Type 7	100380	
	MicroThread	3.5	Type 5	100376	
		4.0	Type 6	100378	
BEGO IMPLANT SYSTEMS		4.5, 5.0	Type 7	100380	
	Semados	S3.25 - S5.5, RI 3.75 - RI 5.5	Type 26	100425	
	Semados Mini-Implant	2.7, 2.9, 3.1	Type 2	100355	
	BIOHORIZONS	Internal, Tapered Internal, Single-stage	Yellow platform	Type 32	100440
		Green platform	Type 27	100431	
		Blue platform	Type 33	100442	
External		3.5	Type 30	100436	
		4.0-6.0	Type 1	100353	
		Laser-Lok	3.0	Type 26	100425
BIOLOK	External hex	3.45, 6.5	Type 1	100353	
BRAINBASE	MYTIS Arrow S series	B, C, E-Type	Type 47	100478	
		BW, EW-Type	Type 4	100350	
BRENT	blueSKY		Type 49	100480	
BTI IMPLANTS	Externa®	Universal, Universal Plus	Type 1	100353	
	Interna®	Universal, Universal Plus, Wide	Type 27	100431	
	Tiny		Type 9	100368	
CAMLOG		3.3	Type 28	100432	
		3.8, 4.3	Type 23	100411	
		5.0, 6.0	Type 24	100412	
	CHAORUM	External	3.3	Type 1	100353
	Internal	4.8	Type 47	100478	
COCHLEAR BAHA	BI 300		Type 30	100436	
	Flange Fixture Hexagonal connection		Type 9	100368	
COWELLMEDI	Atlas External	3.3	Type 2	100355	
		4.1, 5.0	Type 1	100353	
	Atlas Submerged		Type 21	100402	
DENTAURUM IMPLANTS	tioLogic	3.3-5.5	Type 45	100472	
DENTEGRIS	Straight & Tapered	3.3	Type 26	100425	
	Straight & Tapered	3.75, 4.5, 4.75	Type 32	100440	
	Straight & Tapered	5.5	Type 47	100478	

Implant System	Model	Size	SmartPeg	Article no.
DENTIS DENTAL IMPLANT SOLUTION	Internal		Type 4	100350
	Save Submerged		Type 7	100380
DENTIUM	Implantium	3.4, 3.8, 4.3	Type 6	100378
	Super Line	4.0 - 7.0	Type 7	100380
	Simple Line	3.8, 4.3, 4.8	Type 7	100380
DENTSPLY FRIADENT	Ankylos	A, B, C	Type 16	100388
	Frialit-2	3.4, 6.5	Type 21	100402
	Xive	3.0	Type 22	100404
		3.4 - 5.5	Type 21	100402
DIO	TG		Type 8	100366
	External FSN	3.3	Type 2	100355
		3.75	Type 1	100353
	External FTN	5.0	Type 3	100357
	SM Submerged	3.8	Type 49	100480
		4.5, 5.3	Type 47	100478
DYNA	Dyna Helix DC	3.2, 3.6	Type 27	100431
	Dyna Helix ST	3.6, 4.2, 5.0	Type 27	100431
ECKERMANN LABORATORIUM	All spiral		Type 17	100390
	External-Hex		Type 1	100353
	Triplo		Type 49	100480
EUROTEKNIKA		3.75 - 6.0	Type 1	100353
GC IMPLANT	GENESIO		Type 21	100402
	SETIO		Type 1	100353
IDI EVOLUTION		3.75 - 4.0	Type 18	100398
IMTEC	Hexed-head	3.75 - 4.75	Type 1	100353
INTRA-LOCK	Internal	SQ	Type 5	100376
		Standard, Wide	Type 47	100478
JDENTALCARE	JDÉvolution	3.7, 4.3, 5.0	Type 32	100440
		3.2	Type 38	100455
JMM	POI 2Piece	3.2, 3.7, 4.2	Type 18	100398
	POIEX HACEX	3.7	Type 15	100386
	POIEX HACEX	4.2	Type 49	100480
	POIEX HACEX	4.7, 5.2	Type 47	100478
KEYSTONE DENTAL	Renova	SDI	Type 26	100425
		RDI	Type 27	100431
	Restore	3.3	Type 30	100436
		3.75, 4.0	Type 1	100353
		5.0, 6.0	Type 3	100357
	Prima Connex	SD	Type 27	100431
		RD, WD	Type 33	100442
KLOCKNER IMPLANT SYSTEM	Essential Cone		Type 49	100480
	NK2, S3M, S4, S6, SK2		Type 30	100436
LASAK	Impladent	2.9	Type 2	100355
		3.7	Type 36	100447
		5.0	Type 37	100448
MEDENTIKA	M-Implant	3.5-5.0	Type 38	100455
MEDIGMA TECHNOLOGIES	Fix-A-Dent, FixTite, FixTites		Type 49	100480
	Fix 1		Type 17	100390
MEGAGEN	Rescue External	7.0	Type 3	100357
	EZ Plus Internal	4.0, 4.5, 5.0	Type 7	100380
	EZ Plus Internal	Small	Type 38	100455
	EZ Plus Internal	Re/Wi	Type 7	100380
MICRODENT	MS, MST	3.5	Type 30	100436
	MK	5.6	Type 3	100357
	Trylogic	3.8, 4.3	Type 27	100431
	Universal	4.1, 5.1	Type 1	100353
MIS	Internal Hexagon		Type 32	100440
	Internal Octagon		Type 4	100350
MOZO-GRAU	MG Osseous	Mini, Standard	Type 1	100353
		Maxi	Type 3	100357
	MG Inhex	Standard	Type 38	100455
		Maxi	Type 7	100380

Implant System	Model	Size	SmartPeg	Article no.
NEMRIS	Aesthura Classig	3.75, 4.75	Type 35	100446
NEOBIOTECH	Neoplant External	3.75, 4.0, 5.0, 6.0	Type 1	100353
NEODENT	Alvim CM	3.5-5.0	Type 16	100388
	Alvim Ti	4.3	Type 1	100353
	Alvim II Plus	3.75-5.0	Type 21	100402
	Titamax CM	3.5-5.0	Type 16	100388
	Titamax GT	3.7-4.5	Type 9	100368
	Titamax Ti	3.75, 4.0	Type 1	100353
	Titamax Ti	4.5, 5.0	Type 3	100357
	Titamax II Plus	3.75-5.0	Type 21	100402
	Neoss Implant	3.5-5.5	Type 18	100398
		3.25	Type 38	100455
NOBEL BIOCARE	Active	NP	Type 21	100402
		RP	Type 47	100478
NOBEL DIRECT	Brånemark system, Nobel Speedy Groovy	NP	Type 2	100355
		RP	Type 1	100353
		WP	Type 3	100357
		NP	Type 8	100366
		RP	Type 9	100368
		WP	Type 9	100368
		NP	Type 10	100370
		RP	Type 11	100372
		WP	Type 14	100374
		NP	Type 12	100359
		RP	Type 13	100361
		WP	Type 13	100361
		6.0	Type 40	100462
		RP	Type 1	100353
	WP	Type 3	100357	
OSSEOLINK		3.5-5.0	Type 21	100402
OSPOL IMPLANT SYSTEM		NP, RP	Type 21	100402
	USII, USIII	Mini	Type 2	100355
OSSTEM IMPLANTS		Regular	Type 1	100353
		Wide, Ultra Wide	Type 3	100357
		USIV Conical Mini	Type 8	100366
		GSII Mini	Type 5	100376
		Regular, Wide, Ultra Wide	Type 6	100378
		GSIII Mini	Type 5	100376
		Regular	Type 6	100378
		SSI, SSII, SSIII Regular	Type 4	100350
		Wide	Type 4	100350
			Type 10	100370
OSTEOPLANT	Two piece		Type 10	100370
OSTEO-TI		3.25 - 5.5	Type 26	100425
P-I BRÅNEMARK PHILOSOPHY	Internal Grip Amplified	4.3, 5.1	Type 51	100482
PLAN 1 HEALTH	Hexcel IHC, IHT, Internal Hex		Type 26	100425
	Hexcel EHC, External Hex		Type 1	100353
PLATON	Type I, Type II		Type 4	100350
	Type IV Pro		Type 7	100380
PRODENT ITALIA	Perform SM	3.6, 4.0	Type 49	100480
		4.5, 5.0	Type 7	100380
REVOIS		3.8, 4.3, 5.0	Type 26	100425
SIC IMPLANT SYSTEM	SIC ace	3.4, 4.0, 5.0	Type 44	100466
	SIC pro	3.4, 4.0, 5.0	Type 44	100466
S.I.N. IMPLANT SYSTEM	External Hex	3.4, 4.1, 5.0	Type 1	100353
	Internal Hex	3.8, 4.5, 5.5	Type 15	100386
SOUTHERN IMPLANTS	External Hex		Type 1	100353
	IT-Connection		Type 4	100350
	Tri-nex	3.5	Type 12	100359
	Tri-nex	4.3, 5.0	Type 13	100361
	Extra-Oral		Type 1	100353
STERNGOLD DENTAL	Stern IC Internal Octagon	Regular Neck 4.1	Type 4	100350
	External Hexagon	3.75, 4.0, 5.0 RP	Type 1	100353